



USER MANUAL AT-599I

---

Thank you for purchasing the AT-599I Massage Chair.  
Before using your massage chair, please read the contents of this  
instruction manual and refer to it when needed.

## Table of Contents

Luxury Features .....	2
Safety Precautions .....	3 - 4
Chair Components .....	5 - 6
Chair Installing .....	7 - 9
Moving the Chair .....	10
Grounding Instructions .....	11
Chair Specifications .....	12
Operation Guide .....	13-17
Troubleshooting .....	18
Warnings .....	18

---

### **AT-599I Luxury Massage Chair**

*The Luxury Massage Chair is controlled by a microcomputer that performs a variety of intelligent massage techniques and combinations that are designed to soothe sore muscles, alleviate stress and fatigue, increase metabolism and energy, and improve blood circulation.*

## **01 Luxury Features**

- Auto-scanning automatically or manually customizes a personal massage specifically for each user.
- Invigorating air pressure system in seat and calves.
- Six pre-programmed massage sessions with customizable speed, width and intensity.
- Ability to manually program any combination of massage techniques and to adjust rollers to desired location for kneading and Tapping techniques.
- Fifteen-minute program duration with automatic shut off.
- User-friendly LCD remote controller.
- Ergonomically designed for optimum comfort.

---

## 02 Safety Precautions



When using the massage chair, basic precautions should always be followed to reduce the risk of electric shock, burns, fire or injury to persons. Please read the following precautions before using the massage chair.

- Make sure the power plug is securely inserted into a proper electrical outlet to reduce the risk of short-circuit and fire.
- Always unplug the massage chair from the electrical outlet immediately after using and before cleaning.
- DO NOT use the massage chair in areas where the humidity is high such as a sauna or steam room.
- DO NOT use the massage chair outdoors and avoid areas with direct sunlight.
- If the electricity is off, turn the power switch to the OFF position and unplug the massage chair from the electrical outlet to prevent any damages to your chair when the electricity comes back on.
- Never plug/unplug the power cord plug socket with wet hands.
- Close and continuous supervision is necessary when the massage chair is used by, on or near children or disabled persons.
- No more than one person should use the massage chair at a time.
- The massage chair is intended to be used as described in the Instruction Manual. Do not use any attachments other than those recommended.
- If the supply cord is damaged, it must be replaced by a special cord or assembly available from the manufacturer or its service agent.
- DO NOT attempt to stand or sit on the backrest or footrest in order to prevent accident, injury or damage to the chair.
- Avoid falling asleep while using the massage chair.
- Always keep the massage chair clean, and never insert any objects into or around the gaps of the massage chair.
- DO NOT move the massage chair by pulling on the power supply cord.
- Make sure there are no children, pets or other obstacles behind or underneath the massage chair.
- If you find any rips or tears in the fabric of the massage chair contact our service department.

- 
- This massage chair is designed for household use and should not be used in non-residential settings.
  - Please consult your doctor if you are in doubt as to whether or not you should use this massage chair.
  - Please consult a doctor before using the massage chair if you are pregnant, undergoing any medical treatment, using a medical electronic device such as a pacemaker, or are experiencing any of the following: Malignant tumor, heart disease, back pain, abnormal or curved back bone, osteoporosis or acute illness.
  - If you experience pain while using the massage chair, stop using it immediately and consult your doctor.
  - Never use a massage function on swollen or inflamed parts of the body.
  - DO NOT use the massage chair for more than 30 minutes at a time.
  - DO NOT use the massage chair where aerosol (spray) products are being used or where oxygen is being administered.
  - This massage chair is a non-professional product designed to provide a comfortable massage and should not be substituted for appropriate medical treatment.


---

## 03 Chair Components



---

## MASSAGING POSITION




---

## 04 Chair Installation

### Step 1:


When setting up the chair, make sure there is at least two feet of space behind the backrest and two feet of space in front of the footrest.



### Step 2:

Insert the wire plugs into the seat as demonstrated in PIC 1.

Note: Insert black to black and gray to gray, otherwise damage to the chair may result.






---


**Step 3:**

Unfold the backrest until you hear a “click” sound.



**Step 4:**



Attach the backrest pad and the pillow pad with the zipper, hook and loop fasteners.



---

**Step 5: (Optional)**

Remove the spare pad located in the backrest cover for a deeper massage.




---

## 05 Moving the Chair

1. Restore the footrest to its original position.
2. Hold down on the top of the backrest and tilt approximately 45 degrees. Move the chair using the wheels as the fulcrum.

Note:

- 1) Always return the footrest to its original length and restore it to its original position before moving the chair. Otherwise, injury during movement may result.
- 2) Do not move the chair while a person is seated on it.
- 3) Place a mat on the floor and gently roll the chair over it to avoid damaging the floor.



---

## **06** Grounding Instructions

Fully insert the plug into an appropriate outlet.

An improper equipment-grounding conductor may result in electric shock. Contact a qualified electrician if you are not sure that the product is properly grounded. Do not modify the plug that is included with the massage chair. If the plug does not fit the outlet, contact a qualified electrician to install an outlet to fit the massage chair plug.

---

## 07 Chair Specifications

Description	Specifications
Model No.	AT-599I
Rated Voltage	120V AC
Rated Frequency	60Hz
Rated Power Consumption	180W
Timer	15 minutes
Max Time Rate	60 minutes
Dimensions(LxWxH)	Shipping dimensions: 49.3" x 31" x 34.3"
	Reclined dimensions: 65" x 28.7" x 33.5"
	Upright dimensions: 45" x 28.7" x 49"
Weight	Gross weight: 207 lbs.
	Net weight: 163 lbs.
Length of Wire	Controller wire: 1.3m
	Power supply wire: 2.2m
Usage Condition	Temperature: 10°~ 40°
	Humidity: 30-85 RH
Storage Condition	Temperature: 20°~ 60°
	Humidity: 30-85 RH

## 08 Operation Guide



---

## Massage Chair Control Buttons

### 1. ON/OFF

Massage chair power switch

Press this key once will turn on the backlight of the LCD screen. Press any button to begin the massage functions.

Press this key again, the massage chair stops operating and resets automatically to enter into a standing mode.

### 2. Pause

When the chair is in massage mode and the user wants to temporarily pause it, press the Pause button and all massage functions will stop. Press it again to immediately cease the "Pause" function and resume the original massage function.

### 3/4. Recline/Rise

Massage Angle Adjustment Control Buttons

When plugging in, the massage chair enters into ready state. At the moment you can adjust the massage chair angle freely by pressing the RASE/RECLINE key to enjoy more ideal massage, no matter whether the massage machine is running or not. Press button once to adjust the backrest and footrest to a comfort position.

Press the button again to confirm the position.

### 5/6. Recline/Rise

When plugging in, the massage chair enters into ready state. At the moment you can adjust the footrest angle freely by pressing the RASE/RECLINE key to enjoy more ideal massage, no matter whether the massage machine is running or not.

Press button once to adjust the footrest to a comfort position.

Press the button again to confirm the position.

These two keys are used to adjust footrest angle with intelligent electromotion for achieving more ideal and comfortable massage angle.

The footrest will move from 5 to 90 degree when pressing RASE/RECLINE key until to your ideal angle and then undo.

The massage chair will keep all massage function when doing the adjustment.

---

## **AUTO PROGRAM**

Six pre-programmed massage session are activated by pressing the following buttons: REJUVENATE, REFRESH, DESTRESS, SMART, SHIATSU, AUTOAIR.

### **Manual Back Massage Programming**

#### **7. KNEADING & TAPPING**

Kneading & Tapping Start/Stop Control Button

Press this key once to activate the kneading & tapping function for symmetrical and vigorous massage action on the body's reflex points.

Press this key again to exit this function.

Press the FAST/SLOW button to adjust the speed of massage while the Shiatsu function is activated.

#### **8. KNEADING**

Kneading Function Start/Stop Control Button

Press this key once to activate the kneading function for vigorous kneading action on the body's reflex points.

Press this key again to exit this function.

Press the FAST/SLOW button to adjust the speed of massage while the kneading function is activated.

#### **9. TAPPING**

Tapping Function Start/Stop Control Button

Press this key once to activate the tapping function for vigorous tapping action on the body's reflex points.

Press this key again to exit this function.

Press the WIDTH button and Speed button to adjust the TAPPING width and speed of massage while the tapping function is activated.

#### **10. ROLLING**

Rolling Function Start/Stop Control Button

Press this key once to activate the rolling function for vigorous rolling action on the body's reflex points.

Press this key again to exit this function.

Press the WIDTH button to adjust the width of massage while the Rolling function is activated.



---

## **11/12. FAST/SLOW**

### Speed Adjustment Buttons

Press these two keys to activate the intensity adjustment function on the relative massage action to control the massage intensity of current mode.

Press these two keys to switch from fast to slow or slow to fast with circular adjustment.

## **13/14/15. WIDTH**

### Width Adjustment Buttons

Under the mode of Tapping and Rolling, press this key once to switch massage width within Wide/Medium/Narrow to satisfy different users' demands. But it is invalid under Kneading mode.

## **16/17. UP/DOWN**

### Position Adjustment Button for Spot or Partial Massage

These keys are used for adjusting backrest massage position to achieve spot massage.

Press UP and DOWN keys continuously until the massage wheel rolls to ideal position and then undo the key to get the spot massage function.

Press ROLLING key to eliminate spot function.

## **18. PARTIAL**

### Partial Back Massage

Backrest massage Position Selection key is your convenient choice for backrest spot massage.

Press this key to massage within  $\pm 8$ cm up-and-down under current massage mode.

Press this key again to exit this function.

## **19. SPOT**

### Spot Back Massage

The backrest spot massage choice key is your convenient choice for backrest spot massage.

Press this key to massage under current massage mode.

Press this key again to exit this function.

---

## **20. VIBRATE (SEAT)**

Seat Vibration Control Button

Press this key once to activate seat vibrating massage function for comfortable vibrating massage enjoyment.

Press this key again to exit the function.

## **21. KNEAD (CALVES)**

Calves kneading Function Control Button

Press this key once to activate calves kneading massage function for comfortable kneading massage enjoyment.

Press this key again to exit the function.

## **Setting Air Squeeze**

### **22. AUTO 1**

Air Squeeze Control Button— Auto Program 1

Press this key to activate the seat and calves auto air squeezing program No.1 for comfortable waving massage. Press this key again to exit the function.

### **23. AUTO 2**

Air Squeeze Control Button— Auto Program 2

Press this key to activate the seat and calves auto air squeezing program No.2 for comfortable squeezing massage. Press this key again to exit the function.

### **24. INTENSITY(Seat and Calved)**

Press this key once to increase the intensity of seat & calves air massage

Press again to circulate from weak to strong.

Press this key once to decrease the intensity of seat & calves air massage

Press again to circulate from strong to weak.

---

## 09 Troubleshooting & Product Maintenance

This massage chair is controlled by a microcomputer and is susceptible to power surges of +10 to -10%.

The following situations may display in the event of a power surge:

Situation	Solution
The massage chair is connected to the power supply but is not functioning.	<ol style="list-style-type: none"><li>1. Check to see that the plug is well connected and working.</li><li>2. Make sure that the massage chair switch is turned to the "I" position.</li></ol>
The remote control works, but the massage chair is not functioning.	The massage chair may have overheated due to prolonged use. Turn off the power and allow the chair to cool for 30-50 minutes.
Both the remote control and the massage chair are not functioning.	<ol style="list-style-type: none"><li>1. Make sure the chair is plugged into a working outlet.</li><li>2. Turn the power OFF and then ON again.</li></ol>

## 10 Warnings

- \* Do not use water and other cleaning detergents on the massage chair.
- \* Do not stack, stab, slit or harshly pull on the surface of the massage chair.
- \* Gently press the buttons during operation and avoid hitting the buttons with hard objects.
- \* Pull out the power supply plug after each use.
- \* Do not use the massage chair for more than 30 minutes at a time.

